

# SO-850

## Slit Lamp Module

...optics for life!

**Robust**

**Lightweight**

**Precision engineered**

**Left/right selectable eyepiece**

**Selectable cobalt blue filter**

The SO-850 Slit Lamp Module is a portable slit lamp compatible with Welch Allyn and Heine TL battery handles.

It is suitable for high-resolution examination of corneal surfaces, the corneal stroma and the anterior chamber.

The SO-850 is the perfect addition to any Welch Allyn or Heine TL diagnostic set. The eyepiece may be swivelled for comfortable observation with either eye. The eyepiece may be focussed to compensate for refractive error.

A cobalt blue filter may be selected for enhanced viewing where fluorescein is used.



### SPECIFICATIONS: SO-850 Series

**Viewing System:**

Monocular

**Working Distance:**

40mm

**Magnification:**

6 X

**Filter:**

Selectable cobalt blue filter

**Lamp:**

WA 4900 3.5V 0.8A

**Slit Angle:**

35° left or right

**Slit:**

0.4 x 12mm



E-mail: [office@eye4vision.com](mailto:office@eye4vision.com)  
Web site: [www.eye4vision.com](http://www.eye4vision.com)

