



0481

KERATACx™ Series ICRs & Vejarano Delaminator System™




ICR SERIES :	KERATACx™	KERATACx Plus™	KERATACx SP™	KERATACx V™ *
Model Number :	CR-45 / CR-90 CR-120 / CR-160 CR-210 / CR-320	CRP-45 / CRP-90 CRP-120 / CRP-160 CRP-210 / CRP-320	CRS-45 / CRS-90 CRS-120 / CRS-160 CRS-210 / CRS-320	CRV-45 / CRV-90 CRV-120 / CRV-160 CRV-210 / CRV-320
Optical Zone :	4.8mm - 5.2mm	5.4mm - 5.6mm	5.8mm - 6.0mm	6.0mm - 6.1mm
Peak Style :	Rounded, Off-center	Rounded, Centered	Rounded, Centered	Rounded, Centered
Segment Height :	100µ - 450µ in 0.25µ increments	100µ - 450µ in 0.25µ increments	100µ - 450µ in 0.25µ increments	100µ - 450µ in 0.25µ increments
Segment Width :	0.62mm	0.62mm	0.62mm	0.62mm
Degrees :	45° / 90° / 120° 160° / 210° / 320°	45° / 90° / 120° 160° / 210° / 320°	45° / 90° / 120° 160° / 210° / 320°	45° / 90° / 120° 160° / 210° / 320°
Holes :	2 holes, 1 per end @ 0.23mm	2 holes, 1 per end @ 0.23mm	2 holes, 1 per end @ 0.23mm	2 holes, 1 per end @ 0.23mm
Stromal Depth :	70% - 75%	70% - 75%	70% - 75%	70% - 75%
ICR Material :	Clear Perspex C.Q. UV - PMMA	Clear Perspex C.Q. UV - PMMA	Clear Perspex C.Q. UV - PMMA	Clear Perspex C.Q. UV - PMMA

* KERATACx V™ intra corneal rings are used when 4 segments are required.

Note : An optional Delaminator to insure a proper tunnel radius must be purchased for KERATACx™, KERATACx 'SP'™ & KERATACx 'V'™ . (Recommended)

Vejarano Delaminator System™




The World's First affordable corneal delaminating system for Keratoconus, designed to improve results and prevent corneal perforations, and is manufactured of the finest medical grade Titanium or Stainless steel.

Surgeons around the world report that patients experience NO glare, Halos, Sparkles or Star bursts, after implantation of KERATACx™ branded segments with the Vejarano Delaminator System™

Presented in its own 'State of the Art' sterilization tray

the Vejarano Delaminator System™ includes:

- Visual Axis Marker - 1 Unit
- Corneal Marker - 1 Unit
- Pre Delaminator - 1 Unit
- 'Y' Spatula - 1 Unit
- Delaminator - 2 Units (Left & Right)
- Inserter - 2 Units (Left & Right)
- Extractor - 2 Units (Left & Right)
- Forceps - 1 Unit
- 'V' Glide Automatic Vacuum System - 1 Unit



Vejarano Delaminator - Platinum Upgrade Set for KERATACx™ Intra Corneal Segments

Convert your Basic Set (KERATACx Plus™) with the Upgrade Set for KERATACx™ 'CR' ICRs.

The Upgrade Set contains :

- 5mm Corneal Marker - 1 unit
- Delaminator for KERATACx™ 'CR' Rings - 2 Units (Left & Right)
- 'Y' Spatula - 1 unit



Note:

Delaminators for KERATACx™, KERATACx 'SP'™ & KERATACx 'V'™, Instruments and Markers can be purchased separately. Please email us for pricing.

Warranty

The Vejarano Delaminating System™ is supplied in an auto-clave ready tray, with easy to follow USFDA, CE and ISO approved cleaning instructions. Each system is guaranteed against manufacturing defects and workmanship. If the Delaminator becomes bent or misaligned, repairs are free of charge for the first repair, exclusive of all shipping and handling.



Contact us for more information

Email : office@eye4vision.com

Web : www.eye4vision.com