

SO - 5900

S P E C I F I C A T I O N S

OPTICAL HEAD

Viewing System	Binocular stereoscopic, convergent angle 10° for viewing comfort. Image erection by high resolution prism system.
Eyepieces	Tube inclination 45° High field number 22 for wide field of view High eye relief 24.6mm and reversible eye cups for ease of use by spectacle wearers.
Magnification	Zoom, 4.2x to 25x continuous.
Field of View	53mm to 9mm depending on magnification.
Refractive Error Correction	±5D anisometropia
Pupillary Distance	50mm to 76mm, eyepiece tubes locked for ease of adjustment.
Focusing	Fine focus range ±25mm.
Mould Protection	Anti-fungus pellets (three year protection) standard.

ILLUMINATION

Type	Coaxial, direct optical delivery system.
Lamp	20W LED, 65,000 lux, lifetime of the microscope.
UV and Blue light hazard	UV and blue light hazard: Filtered to 435nm
Colour Temperature	4,100°K
Intensity range	10 levels

XY

Number of axis	Two
Range	± 25mm from centre
Return to centre	On start up or via a no touch sterile, infra red sensor.

FOOT CONTROL

Functions	Focus, Zoom, Light Intensity, XY
Dimension	310 x 210 x 70 mm
Weight	5.1Kg
Cable length	2.0m
Mounting	Yes, to the main pillar

MOUNTING / MECHANICAL

Arm	Elbow and angle poise movement in horizontal and vertical directions. Variable friction locks.
Reach	1150 mm
Head Tilt	+5° - 45°
Vertical Tension	Adjustable gas springs to set lifting force.
Height	1600 mm to top of pillar
Head Rotation	330° about the centre of the XY

POWER SUPPLY

Mains Input	100-240VAC, 50-60Hz, 2.0A max. @ 100VAC
Mains Output Regulation	< ±5%
Protection	Over voltage, over current and short circuit protection
Battery / Auxiliary Input	12V DC, polarity and current surge protected.
Approvals	Medical grade approvals CSA C22.2 No. 60601.1-M90 UL 60601-1, 2nd edition IEC 60601-1 CE marked EN 60601-1-2:2001

PACKAGING AND WEIGHTS

Microscope	810 x 520 x 31 mm 26 kg
Floor stand (two cartons)	900 x 270 x 290 mm 25 kg
Video System (optional)	420 x 490 x 340 mm 5 kg
Nett Weight	SO-5900: 72 kg SO-5900AV: 76 kg

OPTIONAL UPGRADES

Video System	SO-5900V. Includes beam splitter, c-mount adapter, CCD camera, LCD monitor and monitor arms, cables and power supplies
Assistant Microscope	Binocular, non-coaxial assistant microscope for left or right mounting.

SO-5900 - STANDARD
 SO-5900A + ASSISTANT MICROSCOPE
 SO-5900V + VIDEO SYSTEM
 SO-5900AV + ASSIST. MICROSCOPE + VIDEO SYSTEM



Lic. No. QEC 10722

Quality Endorsed Company
 Certified to ISO 9001