

**NYLON - POLYAMIDE 6/6 NON ABSORBABLE SYNTHETIC SUTURE**



Non-absorbable  
 Monofilament

**Spatula / Lancet tip**



**CE 1275**  
 Made in Germany




Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)



A007PA	<b>11/0</b>	15	3/8	4.2 100	S
A090PA		20	1/2	5.0 150	
A088PA		20 / 30	1/2	5.6 150	S

A008PA	<b>10/0</b>	15	3/8	4.2 150	S
A073PA		45		4.8 150	
A052PA		30		4.8 150	S
A009PA		30		4.8 150	S
A058PA		45		4.8 150	S
A010PA		30		4.8 150	
A047PA		38		4.8 150	
A067PA		38		4.8 200	
A055PA		45		4.8 200	
A075PA		15		5.6 150	S
A061PA		30		5.6 150	S
A051PA		30		5.6 150	
A037PA		15		6.0 100	
A039PA		30		6.0 100	
A092PA		15		6.0 150	S
A062PA		30		6.0 150	
A066PA		38		6.0 150	
A102PA		30		6.0 200	
A024PA		4		6.2 150	S
A043PA		15		6.2 150	S
A013PA		20		6.2 150	S
A014PA		30	3/8	6.2 150	S
A069PA		38		6.2 150	S
A087PA		10		6.2 150	
A085PA		15		6.2 150	
A034PA		20		6.2 150	
A016PA		30		6.2 150	
A036PA		38		6.2 150	
A045PA		45		6.2 150	
A048PA		75		6.2 150	
A091PA		15		6.4 150	S
A080PA		30		6.4 150	
A065PA		45		6.4 150	
A076PA		15		6.4 200	S
A025PA		30		6.4 200	S
A020PA		38		6.4 200	S
A074PA		45		6.4 200	S
A026PA		30		6.4 200	
A023PA		38		6.4 200	
A041PA		45		6.4 200	
A053PA		75		6.4 200	
A057PA		38		6.4 220	
A060PA		38		6.6 200	
A070PA		30	3/8	7.0 200	

3/8 Circle



1/2 Circle



**NYLON - POLYAMIDE 6/6 NON ABSORBABLE SYNTHETIC SUTURE**



Non-absorbable

Monofilament

**Spatula / Lancet tip**



**CE 1275**

Made in Germany



Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)

A054PA	<b>10/0</b>	13	1/2	5.6 100	S
A098PA		15		5.2 150	
A040PA		15		5.6 150	S
A001PA		20		5.6 150	S
A002PA		30		5.6 150	S
A004PA		30	1/2	5.6 150	
A071PA		30		6.0 150	
A101PA		30		6.2 150	
A081PA		23		6.4 200	S
A068PA		45		6.4 200	
A096PA		30		7.1 150	
A077PA		30	1/2	8.0 200	

A029PA	<b>10/0</b>	15	50-90°	4.8 150	S
A031PA		20		4.8 150	
A032PA		30		4.8 150	
A089PA		30		5.0 150	
A104PA		38	50-90°	6.4 200	
A097PA		30	Straight	16.0 150	S

A033PA	<b>9/0</b>	20	50-90°	4.8 150	
A027PA		30	1/4	6.4 200	
A083PA		45	1/4	6.5 200	
A030PA		15	3/8	4.8 150	S
A011PA		30		4.8 150	
A038PA		15		6.0 100	
A072PA		30		6.0 150	
A015PA		20		6.2 150	S
A079PA		30		6.2 150	S
A017PA		30		6.2 150	
A093PA		30		6.4 150	
A018PA		30		6.4 200	S
A044PA		38		6.4 200	S
A021PA		30		6.4 200	
A042PA		38		6.4 200	
A049PA		45		6.4 200	
A050PA		75		6.4 200	
A086PA		75	3/8	8.0 200	
A094PA		20	1/2	4.8 150	
A003PA		20		5.6 150	S
A005PA		20		5.6 150	
A056PA		30		5.6 150	
A078PA		30	1/2	8.0 200	

1/4 Circle



3/8 Circle




1/2 Circle



Bi-curved 50-90°



Straight



**NYLON - POLYAMIDE 6/6 NON ABSORBABLE SYNTHETIC SUTURE**



Non-absorbable

Monofilament

**Spatula / Lancet tip**



**CE 1275**

Made in Germany



Code	Suture		Needle		
	USP	Length ( cm )	Curve	Length / Diam. ( mm / micron )	Armed ( S-Single )

A012PA	<b>8/0</b>	30	3/8	4.8 150	
A006PA		30	1/2	5.6 150	
A084PA		45	1/4	6.0 200	
A028PA		30	1/4	6.4 200	
A019PA		30	3/8	6.4 200	S
A035PA		38		6.4 200	S
A059PA		45		6.4 200	S
A022PA		30		6.4 200	
A046PA		38	3/8	6.4 200	

A095PA	<b>7/0</b>	30	3/8	6.4 200	
A082PA		23	1/2	10.0 230	S

A063PA	<b>6/0</b>	45	1/4	8.0 300	
--------	------------	----	-----	---------	--

A064PA	<b>5/0</b>	45	1/4	8.5 350	
A103PA		45	3/8	10.0 350	

A099PA	<b>4/0</b>	45	1/4	8.5 450	
A100PA		75	3/8	8.5 450	S



1/4 Circle



3/8 Circle



1/2 Circle



**NYLON - POLYAMIDE 6/6 NON ABSORBABLE SYNTHETIC SUTURE**



Non-absorbable

Monofilament

**Reverse Cutting needle**



**CE 1275**

Made in Germany



Code	Suture		Needle		
	USP	Length ( cm )	Curve	Length / Diam. ( mm / micron )	Armed ( S-Single )

D006PA	<b>6/0</b>	45	3/8	7.8 200	S
D021PA		45		16.0 430	S
D024PA		45		19.0 430	S
D022PA	<b>5/0</b>	45		16.0 430	S
D025PA		45		19.0 430	S
D023PA	<b>4/0</b>	45		16.0 430	S
D026PA		45	3/8	19.0 430	S



This is an extract from our delivery program, and we also consider any individual requests for needle type, suture material, suture lengths and needle-suture combinations, as well as packaging requests.



E-mail : [office@eye4vision.com](mailto:office@eye4vision.com)

Website : <http://www.eye4vision.com/>