



















General Closure

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 016 PA	10/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
A 051 PGA	8/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	3/8	single
A 030 PGA	8/0	PGA	45 cm	Spatula		0,20 mm	6,4 mm	3/8	double
E 007 PGA	8/0	PGA	30 cm	Taperpoint		0,15 mm	6,4 mm	3/8	single
A 065 PGA	7/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	3/8	single
A 034 PGA	7/0	PGA	45 cm	Spatula		0,20 mm	6,4 mm	3/8	double
A 064 PGA	7/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	1/2	double





Glaucoma

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 016 PA	10/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
A 017 PA	9/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
D 112 SS	6/0	Silk black	30 cm	Reverse Cutting		0,23 mm	8,0 mm	3/8	single












Cataract

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 043 PA	10/0	Nylon black	15 cm	Spatula		0,15 mm	6,2 mm	3/8	single
A 016 PA	10/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
A 076 PA	10/0	Nylon black	15 cm	Spatula		0,20 mm	6,4 mm	3/8	single
A 026 PA	10/0	Nylon black	30 cm	Spatula		0,20 mm	6,4 mm	3/8	double
A 026 SB	8/0	Silk blue	45 cm	Spatula		0,20 mm	6,4 mm	3/8	double







Specialities

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 033 PP	10/0	Polypropylene blue	10 cm	Spatula		0,10 mm	16,0 mm 4,0 mm	Straight 3/8	double
A 028 PP	10/0	Polypropylene blue	23 cm	Spatula		0,15 mm	16,0 mm	Straight	double
A 029 PP	10/0	Polypropylene blue	20 cm	Spatula		0,15 mm	5,6 mm	1/2	double

Strabism

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 002 PEG	6/0	Polyester green	45 cm	Spatula		0,33 mm	8,0 mm	1/4	double
A 001 PEG	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
A 008 PEG	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/2	double
A 066 PGA	10/0	PGA	30 cm	Spatula		0,15 mm	6,2 mm	1/2	single
A 067 PGA	9/0	PGA	30 cm	Spatula		0,15 mm	6,2 mm	1/2	single
A 063 PGA	9/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	3/8	single
A 006 PGA	6/0	PGA	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
A 019 PGA	5/0	PGA	45 cm	Spatula		0,34 mm	8,0 mm	1/4	double
A 068 PGA	9/0	PGA	30 cm	Spatula		0,15 mm	6,2 mm	1/2	double
A 069 PGA	6/0	PGA	30 cm	Spatula		0,22 mm	6,4 mm	1/2	double













Retina

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 002 PEG	6/0	Polyester green	45 cm	Spatula		0,33 mm	8,0 mm	1/4	double
A 001 PEG	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
A 008 PEG	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/2	double
A 006 PGA	6/0	PGA	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
A 019 PGA	5/0	PGA	45 cm	Spatula		0,34 mm	8,0 mm	1/4	double

Keratoplasty

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
A 088 PA	11/0	Nylon black	30cm	Spatula		0,15 mm	5,6 mm	1/2	double
A 032 PA	10/0	Nylon black	30cm	Spatula		0,15 mm	4,8 mm	50-90°	double
A 071 PA	10/0	Nylon black	30cm	Spatula		0,15 mm	6,0 mm	1/2	double

Lid Surgery

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
D 033 PP	6/0	Polypropylene blue	45 cm	Reverse Cutting		0,25 mm	12,0 mm	3/8	single
D 034 PP	6/0	Polypropylene blue	45 cm	Reverse Cutting		0,25 mm	12,0 mm	3/8	double
D112 SS	6/0	Silk black	30 cm	Reverse Cutting		0,23 mm	8,0 mm	3/8	single
D 109 SS	5/0	Silk black	45 cm	Reverse Cutting		0,43 mm	19,0 mm	3/8	single
D 054 SS	4/0	Silk black	45 cm	Reverse Cutting		0,43 mm	19,0 mm	3/8	single
D 057 SS	4/0	Silk black	45 cm	Reverse Cutting		0,43 mm	13,0 mm	3/8	single
D 006 PA	6/0	Nylon black	45 cm	Reverse Cutting		0,20 mm	7,8 mm	3/8	single
D 009 PV	5/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	12,0 mm	1/2	single
D 010 PV	4/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	16,0 mm	1/2	single
D 011 PV	5/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	16,0 mm	3/8	single
D 012 PV	4/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	16,0 mm	3/8	single

Ⓢ PGA corresponds to ETHICON "VICRYL", PVDF – an alternative to POLYPROPYLENE

