



















## General Closure

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 016 PA</b>	10/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
<b>A 051 PGA</b>	8/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	3/8	single
<b>A 030 PGA</b>	8/0	PGA	45 cm	Spatula		0,20 mm	6,4 mm	3/8	double
<b>E 007 PGA</b>	8/0	PGA	30 cm	Taperpoint		0,15 mm	6,4 mm	3/8	single
<b>A 065 PGA</b>	7/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	3/8	single
<b>A 034 PGA</b>	7/0	PGA	45 cm	Spatula		0,20 mm	6,4 mm	3/8	double
<b>A 064 PGA</b>	7/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	1/2	double





## Glaucoma

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 016 PA</b>	10/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
<b>A 017 PA</b>	9/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
<b>D 112 SS</b>	6/0	Silk black	30 cm	Reverse Cutting		0,23 mm	8,0 mm	3/8	single












## Cataract

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 043 PA</b>	10/0	Nylon black	15 cm	Spatula		0,15 mm	6,2 mm	3/8	single
<b>A 016 PA</b>	10/0	Nylon black	30 cm	Spatula		0,15 mm	6,2 mm	3/8	double
<b>A 076 PA</b>	10/0	Nylon black	15 cm	Spatula		0,20 mm	6,4 mm	3/8	single
<b>A 026 PA</b>	10/0	Nylon black	30 cm	Spatula		0,20 mm	6,4 mm	3/8	double
<b>A 026 SB</b>	8/0	Silk blue	45 cm	Spatula		0,20 mm	6,4 mm	3/8	double







## Specialities

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 033 PP</b>	10/0	Polypropylene blue	10 cm	Spatula		0,10 mm	16,0 mm 4,0 mm	Straight 3/8	double
<b>A 028 PP</b>	10/0	Polypropylene blue	23 cm	Spatula		0,15 mm	16,0 mm	Straight	double
<b>A 029 PP</b>	10/0	Polypropylene blue	20 cm	Spatula		0,15 mm	5,6 mm	1/2	double

## Strabism

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 002 PEG</b>	6/0	Polyester green	45 cm	Spatula		0,33 mm	8,0 mm	1/4	double
<b>A 001 PEG</b>	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
<b>A 008 PEG</b>	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/2	double
<b>A 066 PGA</b>	10/0	PGA	30 cm	Spatula		0,15 mm	6,2 mm	1/2	single
<b>A 067 PGA</b>	9/0	PGA	30 cm	Spatula		0,15 mm	6,2 mm	1/2	single
<b>A 063 PGA</b>	9/0	PGA	30 cm	Spatula		0,20 mm	6,4 mm	3/8	single
<b>A 006 PGA</b>	6/0	PGA	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
<b>A 019 PGA</b>	5/0	PGA	45 cm	Spatula		0,34 mm	8,0 mm	1/4	double
<b>A 068 PGA</b>	9/0	PGA	30 cm	Spatula		0,15 mm	6,2 mm	1/2	double
<b>A 069 PGA</b>	6/0	PGA	30 cm	Spatula		0,22 mm	6,4 mm	1/2	double













## Retina

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 002 PEG</b>	6/0	Polyester green	45 cm	Spatula		0,33 mm	8,0 mm	1/4	double
<b>A 001 PEG</b>	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
<b>A 008 PEG</b>	5/0	Polyester green	45 cm	Spatula		0,35 mm	8,0 mm	1/2	double
<b>A 006 PGA</b>	6/0	PGA	45 cm	Spatula		0,35 mm	8,0 mm	1/4	double
<b>A 019 PGA</b>	5/0	PGA	45 cm	Spatula		0,34 mm	8,0 mm	1/4	double

## Keratoplasty

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>A 088 PA</b>	11/0	Nylon black	30cm	Spatula		0,15 mm	5,6 mm	1/2	double
<b>A 032 PA</b>	10/0	Nylon black	30cm	Spatula		0,15 mm	4,8 mm	50-90°	double
<b>A 071 PA</b>	10/0	Nylon black	30cm	Spatula		0,15 mm	6,0 mm	1/2	double

## Lid Surgery

Code	Suture			Needle					
	USP	Type	Length	Type		Diameter	Length	Curve	Armed
<b>D 033 PP</b>	6/0	Polypropylene blue	45 cm	Reverse Cutting		0,25 mm	12,0 mm	3/8	single
<b>D 034 PP</b>	6/0	Polypropylene blue	45 cm	Reverse Cutting		0,25 mm	12,0 mm	3/8	double
<b>D112 SS</b>	6/0	Silk black	30 cm	Reverse Cutting		0,23 mm	8,0 mm	3/8	single
<b>D 109 SS</b>	5/0	Silk black	45 cm	Reverse Cutting		0,43 mm	19,0 mm	3/8	single
<b>D 054 SS</b>	4/0	Silk black	45 cm	Reverse Cutting		0,43 mm	19,0 mm	3/8	single
<b>D 057 SS</b>	4/0	Silk black	45 cm	Reverse Cutting		0,43 mm	13,0 mm	3/8	single
<b>D 006 PA</b>	6/0	Nylon black	45 cm	Reverse Cutting		0,20 mm	7,8 mm	3/8	single
<b>D 009 PV</b>	5/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	12,0 mm	1/2	single
<b>D 010 PV</b>	4/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	16,0 mm	1/2	single
<b>D 011 PV</b>	5/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	16,0 mm	3/8	single
<b>D 012 PV</b>	4/0	PVDF blue	45 cm	Reverse Cutting		0,43 mm	16,0 mm	3/8	single

Ⓢ PGA corresponds to ETHICON "VICRYL", PVDF – an alternative to POLYPROPYLENE

