

POLYPROPYLENE - NON ABSORBABLE SYNTHETIC SUTURE



Non-absorbable

Monofilament

Spatula / Lancet tip



CE 1275

Made in Germany



Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)

A033PP	10/0	10	Straight	16.0 100	
			3/8	4.0 100	

A028PP	10/0	23	Straight	16.0 150	
--------	-------------	----	----------	----------	--

A029PP		20	1/2	5.6 150	
--------	--	----	-----	---------	--

A005PP		15	3/8	4.0 100	S
--------	--	----	-----	---------	---

A006PP		30		4.0 150	
--------	--	----	--	---------	--

A007PP		30	3/8	4.0 150	S
--------	--	----	-----	---------	---

A001PP		20	1/2	5.6 150	S
--------	--	----	-----	---------	---

A003PP		30	1/2	5.6 150	
--------	--	----	-----	---------	--

A010PP		30	3/8	6.2 150	S
--------	--	----	-----	---------	---

A011PP		30		6.2 150	
--------	--	----	--	---------	--

A016PP		30	3/8	6.4 200	
--------	--	----	-----	---------	--

A008PP	9/0	30	3/8	4.2 150	
--------	------------	----	-----	---------	--

A002PP		20	1/2	5.6 150	S
--------	--	----	-----	---------	---

A004PP		30	1/2	5.6 150	
--------	--	----	-----	---------	--

A014PP		15	3/8	6.4 200	S
--------	--	----	-----	---------	---

A012PP		30		6.4 200	
--------	--	----	--	---------	--

A017PP		30	3/8	6.4 200	S
--------	--	----	-----	---------	---

A009PP	8/0	30	3/8	4.0 150	S
--------	------------	----	-----	---------	---

A015PP		15		6.4 200	S
--------	--	----	--	---------	---

A013PP		30	3/8	6.4 200	
--------	--	----	-----	---------	--

3/8 Circle




1/2 Circle



Bi-curved 50-90°

(A033PP)

Straight



POLYPROPYLENE - NON ABSORBABLE SYNTHETIC SUTURE



Non-absorbable

Monofilament

Reverse Cutting needle




CE 1275

Made in Germany

Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)

D002PP	5/0	45	3/8	19.0 430	S
--------	------------	----	-----	----------	---

D001PP	4/0	45	3/8	19.0 430	S
--------	------------	----	-----	----------	---



This is an extract from our delivery program, and we also consider any individual requests for needle type, suture material, suture lengths and needle-suture combinations, as well as packaging requests.



E-mail : office@eye4vision.com

Website : http://www.eye4vision.com/