

SILK - NON ABSORBABLE NATURAL SUTURE



Non-absorbable

Braided Multifilament

Spatula / Lancet tip




CE 1275

Made in Germany



Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)

					
A008SB	10/0	15	3/8	4.2 100	
A009SB		30		4.8 150	S
A011SB		30		4.8 150	
A049SB		15		5.6 150	S
A051SB		30		6.0 150	
A052SB		30		6.2 150	S
A044SB		45		6.2 150	
A045SB		45		6.4 200	
A053SB		30		6.4 200	S
A054SB		30	3/8	6.4 200	
A001SB		20	1/2	5.6 150	S
A004SB		30	1/2	5.6 150	

A021SB	9/0	45	1/4	6.4 200	
A012SB		30	3/8	4.8 150	
A050SB		30		6.0 150	
A017SB		30		6.4 200	S
A019SB		30		6.4 200	
A039SB		38		6.4 200	
A041SB		45	3/8	6.4 200	
A002SB		20	1/2	5.6 150	S
A005SB		30		5.6 150	
A038SB		30	1/2	7.1 150	
A022SB		45	1/4	6.4 200	


A022SB	8/0	45	1/4	6.4 200	
A013SB		30	3/8	4.8 150	
A047SB		30		6.0 150	
A065SS		30		6.0 150	
A067SS		45		6.0 150	S
A014SB		30		6.2 150	S
A015SB		30		6.2 150	
A036SB		38		6.2 150	
A010SB		45		6.2 150	
A016SB		45		6.2 150	
A018SB		30		6.4 200	S
A061SB		30		6.4 200	S
A020SB		30		6.4 200	
A029SB		38		6.4 200	S
A032SB		38		6.4 200	
A027SB		45		6.4 200	S
A026SB		45		6.4 200	
A057SS		45		6.4 200	
A042SB		75		6.4 200	
A034SB		75		8.0 200	
A031SB		75	3/8	9.7 230	
A003SB		20	1/2	5.6 150	S
A006SB		30		5.6 150	
A066SS		30		8.0 220	
A060SS		45	1/2	8.0 200	

A046SS	7/0	45	3/8	6.4 200	
A030SB		75		6.4 200	
A059SB		30		8.0 250	
A058SS		75	3/8	8.0 250	
A033SS		45	1/2	8.0 250	

1/4 Circle

3/8 Circle

1/2 Circle



SILK - NON ABSORBABLE NATURAL SUTURE



Non-absorbable

Braided Multifilament

Spatula / Lancet tip



CE 1275

Made in Germany



Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)

A035SS	6/0	45	1/4	6.4 350	
A064SS		45		8.0 340	
A024SS		45	1/4	8.5 350	
A048SS		45	3/8	6.4 200	
A040SS		45		8.0 300	
A055SS		75		8.0 300	
A043SB		75		8.0 340	
A028SS		45		8.5 350	
A037SS		45	1/2	8.0 300	
A007SS		45	1/2	8.5 350	S


A069SS	5/0	20	1/4	8.0 460	S
A025SS		45	1/4	8.5 350	
A056SS		45	3/8	8.0 350	S
A068SS		20	3/8	8.0 460	S

A023SB	4/0	30	1/4	6.0 350	
A063SS		20	1/4	8.0 460	S
A062SS		20	3/8	8.0 460	S

1/4 Circle

3/8 Circle

1/2 Circle



SILK - NON ABSORBABLE NATURAL SUTURE



Non-absorbable

Braided Multifilament

Reverse Cutting needle



CE 1275


Made in Germany

Code	Suture		Needle		
	USP	Length (cm)	Curve	Length / Diam. (mm / micron)	Armed (S-Single)

D009SS	7/0	45	3/8	7.8 200	
D010SS	6/0	45	3/8	7.8 200	
D112SS		30		8.0 230	S
D039SS		75	3/8	16.0 430	S

D040SS	5/0	45	3/8	16.0 430	S
D054SS		45	3/8	19.0 430	S

D057SS	4/0	45	3/8	16.0 430	S
D109SS		45	3/8	19.0 430	S



This is an extract from our delivery program, and we also consider any individual requests for needle type, suture material, suture lengths and needle-suture combinations, as well as packaging requests.



E-mail : office@eye4vision.com

Website : http://www.eye4vision.com/