

CODE	SIZE	DESCRIPTION	PCS. / BOX
D1020	40 x 40cm	Aperture Drape	10
D1021	40 x 40cm	Mini Incise Drape	10
D1035	15 x 20cm	Mini Incise Drape	10
D1061	122 x 127cm	Incise Drape with Pouch	5
D1062	122 x 127cm	Incise Drape with Pouch	5
D1063	122 x 127cm	Aperture Drape with Pouch	5
M0275	25G x 1 1/2"	Standard Retrobulbar Needle	10
M0278	19G x 1"	Sub-Tenon Anaesthesia Cannula	10
M0285	30G x 5/8"	Facial Infiltration Needle	10
M0637	25G x 1 1/2"	Atkinson Point Retrobulbar Needle	10
M0638	23G x 1 1/2"	Atkinson Point Retrobulbar Needle	10
M0641	27G x 7/8"	Peribulbar Needle	10
M0641A NEW	27G x 1 1/4"	Peribulbar Needle	10
M0642	25G x 1 1/4"	Peribulbar Needle	10
M0642A	25G x 7/8"	Peribulbar Needle	10
M0642B	25G x 5/8"	Peribulbar Needle	10
M0643	23G x 1 1/4"	Peribulbar Needle	10
M0644	23G x 7/8"	Peribulbar Needle	10
M2272	26G x 7/8"	Visco-Elastic Cannula (Angled)	10
M2273	30G x 7/8"	Rycroft Anterior Chamber Cannula	10
M2273A	27G x 7/8"	Rycroft Anterior Chamber Cannula	10
M2273B	20G x 7/8"	Rycroft Anterior Chamber Cannula	10
M2273C	27G x 7/8"	Sauter Curved Hydrodissector	10
M2273D	25G x 7/8"	Sauter Curved Hydrodissector	10
M2273E	27G x 7/8"	Helsinki Curved Hydrodissector	10
M2273G	23G x 7/8"	Rycroft Anterior Chamber Cannula	10
M2273H	27G x 7/8"	Curved Hydrodissection Cannula	10
M2275	25G x 7/8"	Rycroft Anterior Chamber Cannula	10
M2639	25G x 1"	Rycroft Anterior Chamber Cannula	10
M2639A	25G x 7/8"	Visco-Elastic Cannula (Angled)	10
M2642	21G x 7/8"	Rycroft Anterior Chamber Cannula	10
M2656	20G x 9/64"	Micro Etched A C Maintainer (Dahan)	5
M2657	20G x 9/64"	Self Retaining A C Maintainer (Lewicky)	5
M2658 NEW	23G x 9/64"	Self Retaining A C Maintainer (Lewicky)	5
M2801	25G x 7/8"	Lin Lasik Irrigating Cannula	10
M2802	25G x 1 1/8"	Buratto Lasik Irrigating Cannula	10
M2804 NEW	25G x 7/8"	Lasik Irrigating Cannula 6 ports	10
M2813 NEW	23G x 7/8"	Lasik Irrigating Cannula	10
M2817	27G x 7/8"	Lasik Irrigating Cannula 4 ports	10
M2818 NEW	27G x 7/8"	Lasik Irrigating Cannula	10
M3000	27G x 5/8"	Capsularhexis Irrigating Cystotome (Straight)	10
M3000D	25G x 5/8"	Capsularhexis Irrigating Cystotome (Formed) (Short curve)	10
M3000E	27G x 5/8"	Capsularhexis Irrigating Cystotome (Formed) (Short curve)	10
M3000F	27G x 5/8"	Capsularhexis Irrigating Cystotome (Formed)	10
M3000G	25G x 5/8"	Capsularhexis Irrigating Cystotome (Formed)	10
M3001	27G x 5/8"	Endocapsular Irrigating Cystotome (Joyce) (Curved)	10

COLOR	GAUGE
Yellow 	30G
Grey 	27G
Brown 	26G
Orange 	25G
Purple 	24G
Blue 	23G
Black 	22G
Green 	21G
Yellow 	20G
Cream 	19G
Pink 	18G
White 	16G

All hub colours comply with ISO6009

KIT PACKERS

All Sterimedix products are available to kit packers, packed in single blisters, non-sterile. Drapes are provided in bulk, non-sterile.

In addition to Sterimedix products we also supply:

- Instrument Wipes
- Eye Shields
- Drains/Wicks

For more information,

Contact us:

contact@eye4vision.com



www.eye4vision.com

CODE	SIZE	DESCRIPTION	PCS. / BOX
M3602	25G x 5/8"	Side Cutting Cystotome (Formed)	10
M3610	25G x 5/8"	Irrigating Cystotome (Formed)	10
M3610A	25G x 5/8"	Irrigating Cystotome (Straight)	10
M3610B	27G x 5/8"	Irrigating Cystotome (Formed)	10
M3610C NEW	25G x 5/8"	Irrigating Cystotome (Formed) (45° Angled Tip)	10
M3610S	25G x 5/8"	Irrigating Cystotome (Formed with short tip)	10
M3612	27G x 5/8"	Irrigating Cystotome (Straight) (Stoke Mandeville)	10
M3615	27G x 7/8"	Irrigating Cystotome (Formed)	10
M3615A	25G x 7/8"	Irrigating Cystotome (Formed)	10
M3626	27G x 5/8"	Irrigating Cystotome (Angled)	10
M3626A	25G x 5/8"	Irrigating Cystotome (Angled)	10
M3627	25G x 5/8"	Irrigating Cystotome (Formed) Short curve	10
M3627A	27G x 5/8"	Irrigating Cystotome (Formed) Short curve	10
M3628B	23G x 5/8"	Irrigating Cystotome (Formed)	10
M3629	25G x 5/8"	Irrigating Cystotome (Formed, Reverse Tip)	10
M3629A	27G x 5/8"	Irrigating Cystotome (Formed, Reverse Tip)	10
M3630	30G x 5/8"	Irrigating Cystotome (Formed)	10
M3899	25G x 7/8"	Visco-Expression Cannula ('U' curved with flat tip)	10
M3899A	25G x 7/8"	Visco-Expression Cannula ('U' curved with round tip)	10
M3900	25G x 7/8"	Hydrodissection Cannula (Angled)	10
M3901	25G x 7/8"	Hydrodissection Cannula ('J' shaped, flattened tip)	10
M3902	27G x 7/8"	Hydrodissection Cannula (Angled)	10
M3904 NEW	27G x 3/4"	Chang Hydrodissection Cannula	10
M3906A	25G x 7/8"	Micro 'J' Shaped Hydrodissector	10
M3907	25G x 7/8"	Tapered Hydrodelineator	10
M3908	25G x 7/8"	Cortical Cleaving Cannula	10
M3909	23G x 7/8"	Visco-Elastic Cannula (Angled)	10
M3910	21G x 7/8"	Visco-Elastic Cannula (Straight)	10
M3910A	21G x 7/8"	Visco-Elastic Cannula (Angled)	10
M3911	27G x 7/8"	Cortical Cleaving Hydrodissector (Angled, Flattened Tip)	10
M3911A NEW	27G x 7/8"	Cortical Cleaving Hydrodissector (Angled)	10
M3926A	27G x 7/8"	Visco-Elastic Cannula (Angled)	10
M3992	30G x 7/8"	Nucleus Hydrodissector	10
M4410	20G	Subretinal Fluid Cannula 6mm (Straight)	10
M4411	20G	Subretinal Fluid Cannula 3mm (Straight)	10
M4442	20G	Self retaining Scleral Set (1x Plug + 1 x Drain)	5
M4443	20G	Self retaining Scleral Set (2 x Plugs)	5
M4450	19G x 1/2"	Silicone Oil Removal Cannula	10
M4451	16G x 3/8"	High Viscosity Injector	10
M4460	20G x 2.5mm	Infusion Cannula x 150mm	10
M4460A	20G x 2.5mm	Infusion Cannula x 300mm	10
M4461	20G x 4.0mm	Infusion Cannula x 150mm	10
M4461A	20G x 4.0mm	Infusion Cannula x 300mm	10
M4462	20G x 6.0mm	Infusion Cannula x 150mm	10

COLOR	GAUGE
Yellow 	30G
Grey 	27G
Brown 	26G
Orange 	25G
Purple 	24G
Blue 	23G
Black 	22G
Green 	21G
Yellow 	20G
Cream 	19G
Pink 	18G
White 	16G

All hub colours comply with ISO6009

KIT PACKERS

All Sterimedix products are available to kit packers, packed in single blisters, non-sterile. Drapes are provided in bulk, non-sterile.

In addition to Sterimedix products we also supply:

- Instrument Wipes
- Eye Shields
- Drains/Wicks

For more information, Contact us:

contact@eye4vision.com



www.eyevision.com

CODE	SIZE	DESCRIPTION	PCS. / BOX
M4462A	20G x 6mm	Infusion Cannula x 300mm	10
M4463	20G x 125mm	High Pressure Extension Tube	10
M4464	20G x 4.3mm	Self Retaining angled Infusion Cannula x 300mm	10
M4465	20G x 4.3mm	Self retaining infusion Cannula x 300mm	10
M4465A	20G x 4.3mm	Self retaining infusion cannula x 150mm	10
M4466	20G x 2.5mm	Self retaining infusion Cannula x 300mm	10
M4466A	20G x 2.5mm	Self retaining infusion Cannula x 150mm	10
M4467	20G x 6.0mm	Self retaining infusion Cannula x 300mm	10
M4467A	20G x 6.0mm	Self retaining infusion Cannula x 150mm	10
M4468	23G x 4.3mm	Infusion Cannula (Straight) x 150mm	10
M4470	20G > 23G	Reducing Cannula	10
M4471	20G x 1 1/4"	Aspiration Cannula	10
M4472	19G x 1 1/4"	Aspiration Cannula	10
M4473	18G x 1 1/4"	Heavy Oil Removal Cannula	10
M4474	20G / 25G	Perfluorocarbon Cannula	10
M4475	20G x 1 1/4"	Silicone Brush	10
M4475A	20G x 1 1/4"	Silicone Brush	10
M4476	25G x 1 1/4"	Biopsy Cannula	10
M4477	25G x 1 1/4"	Aspiration Cannula	10
M4478	23G x 1 1/4"	Aspiration Cannula	10
M4479	19G x 1 1/2"	Flute Handle with Cannula	5
M4480	20G x 1 1/4"	Flute Handle with Cannula	5
M4481	23G x 1 1/4"	Flute Handle with Cannula	5
M4484	20G x 1 1/4"	Active Flute Handle with Cannula	5
M4488	110mm	Backflush Handle	5
M4489	100mm	'Rigid' Backflush Handle	5
M4490	110mm	'Rigid' Backflush Handle + 20G Silicone Brush	5
M4491 NEW	100mm	'Rigid' Backflush Handle + 20G Silicone Brush	5
M4492	100mm	'Rigid' Backflush Handle + 23G Silicone Brush	5
M4493	110mm	Charles Flute Handle	5
M4494 NEW	110mm	Charles Flute Handle + 20G Aspiration Cannula	5
M4495 NEW	110mm	Charles Flute Handle + 23G Aspiration Cannula	5
M4496 NEW	100mm	Charles Flute Handle + 23G Aspiration Cannula	5
M4497 NEW	100mm	'Rigid' Backflush Handle + 20G Aspiration Cannula	5
M4498 NEW	110mm	Heavy Liquid Infusion Cannula	5
M4499 NEW	23G / 30G	Perfluorocarbon Cannula	10
M4595	27G x 1 1/2"	Irrigating Micro Vectis	10
M4619	25G x 1 1/2"	Irrigating Vectis (Single Holed)	10
M4620	25G x 1 1/2"	Irrigating Vectis (Double Holed)	10

COLOR	GAUGE
Yellow	30G
Grey	27G
Brown	26G
Orange	25G
Purple	24G
Blue	23G
Black	22G
Green	21G
Yellow	20G
Cream	19G
Pink	18G
White	16G

All hub colours comply with ISO6009

KIT PACKERS

All Sterimedix products are available to kit packers, packed in single blisters, non-sterile. Drapes are provided in bulk, non-sterile.

In addition to Sterimedix products we also supply:

- Instrument Wipes
- Eye Shields
- Drains/Wicks

For more information, Contact us: contact@eye4vision.com



www.eye4vision.com

CODE	SIZE	DESCRIPTION	PCS. / BOX
M4621	27G x 11/2"	Small Incision Vectis (Single Holed)	10
M4622	25G x 11/2"	Anis Pattern Vectis	10
M5287	30G x 7/8"	Irrigating Iris Hook	10
M5405 NEW	21G x 5/8"	Bimanual Irrigation Cannula	10
M5409 NEW	21G x 5/8"	Bimanual Irrigation Cannula 2 x 0.5mm ports	10
M5414 NEW	21G x 7/8"	Bimanual Irrigation Cannula 2 x 0.5mm ports	10
M5500	25G x 7/8"	Angled Endothelial Spatula	5
M5502	23G x 5/8"	Aspiration Handle 0.3mm side port and micro etched tip (Curved)	5
M5505	21G x 5/8"	Irrigation Handle (Curved)	5
M5509	21G x 5/8"	Irrigation Handle with 2 x 0.5mm ports (Straight) 14 5	5
M5510	20G x 3/4"	Irrigation handle with 2 x 0.5mm ports (Curved)	5
M5512	21G x 5/8"	Aspiration Handle 0.3mm side port, and micro etched tip (Curved)	5
M5514	21G x 5/8"	Irrigation Handle with 2 x 0.5mm ports (Curved)	5
M5515	21G x 7/8"	Aspiration handle with 0.3mm side port (Curved)	5
M5516	21G x 7/8"	Irrigation handle with concave bevelled tip (Curved)	5
M5520	21G x 5/8"	Aspiration Handle thin wall with 0.4mm side port	5
M5521	21G x 5/8"	Irrigating Handle thin wall with 2 x 0.6mm side ports	5
M5522 NEW	21G x 5/8"	Irrigating Handle thin wall with open end	5
M5531 NEW	21G x 5/8"	Bimanual I/A System combining: Irrigating Handle thin wall with 2 x 0.6mm side ports	5 pairs
	21G x 5/8"	Aspiration Handle thin wall with 0.4mm side port	
M5532 NEW	21G x 5/8"	Bimanual I/A System combining: Irrigating Handle thin wall with open end	5 pairs
	21G x 5/8"	Aspiration Handle thin wall with 0.4mm side port	
M5550 NEW	20G	Irrigating Chopper with 2 x 0.6mm ports	5
M5551 NEW	20G	Irrigating Chopper with open end	5
M5552 NEW	20G	Phaco Chopper (Nagahara style)	5
M5580	16G / 21G	I/A Handpiece (Straight tip)	5
M5581	16G / 21G	I/A Handpiece (Curved tip)	5
M5582	16G / 21G	I/A Handpiece (45° angled tip)	5
M5583	16G / 21G	I/A Handpiece (60° angled tip)	5
M5584	16G / 21G	I/A Handpiece (90° angled tip)	5
M5585	16G / 23G	I/A Handpiece (Binkhorst tip)	5
M5620	23G x 1"	Binkhorst Hook (Right)	10
M5621	23G x 1"	Binkhorst Hook (Left)	10
M5622 NEW	25G x 7/8"	Binkhorst Hook (Right)	10
M5623 NEW	25G x 7/8"	Binkhorst Hook (Left)	10
M6524 NEW	21G x 5/8"	Bimanual Aspirating Cannula	10
M6525	21G x 5/8"	Cortex Extractor 0.3mm port	10
M6575	23G x 5/8"	Cortex Extractor 0.3mm port	10

COLOR	GAUGE
Yellow 	30G 
Grey 	27G 
Brown 	26G 
Orange 	25G 
Purple 	24G 
Blue 	23G 
Black 	22G 
Green 	21G 
Yellow 	20G 
Cream 	19G 
Pink 	18G 
White 	16G 

All hub colours comply with ISO6009

KIT PACKERS

All Sterimedix products are available to kit packers, packed in single blisters, non-sterile. Drapes are provided in bulk, non-sterile.

In addition to Sterimedix products we also supply:

- Instrument Wipes
- Eye Shields
- Drains/Wicks

For more information,
Contact us:

contact@eye4vision.com



www.eye4vision.com

CODE	SIZE	DESCRIPTION	PCS. / BOX
M6575A	23G x 5/8"	Cortex Extractor 0.4mm Port	10
M6576	23G x 7/8"	Cortex Extractor 0.3mm port 11 10	10
M6576A	23G x 1/2"	Cortex Extractor 0.3mm port	10
M6594	23G x 5/8"	Cortex Extractor (Charleux)	10
M6595	23G x 1/2"	Cortex Extractor (Charleux)	10
M6596	24G x 1/2"	Cortex Extractor (Charleux)	10
M6600	23G/23G	Standard Simcoe Aspiration through 0.4mm side port	5
M6601	23G/23G	Reverse Simcoe Aspiration through 0.4mm side port	5
M6605	23G/23G	Standard Simcoe Aspiration through 0.3mm side port	5
M6609	23G/23G	Standard Dahan Simcoe (Sandblasted Tip) Aspiration through 0.4mm side port	5
M6611	21G x 1 1/2"	Single Cannula (Straight)	10
M6612	23G x 1 1/2"	Single Cannula (Straight)	10
M6616	21G x 1 1/2"	Twin Cannula (Straight)	5
M6616C	21G x 1 1/2"	Twin Cannula (Curved)	5
M6634	18G/21G	Coaxial I/A Cannula (Straight)	5
M6634B	18G/21G	Coaxial I/A Cannula (Angled)	5
M7500	3.5mm / 4.0mm	Scleral Marker	10
M8603	21G x 7/8"	Capsule Polisher	10
M8604	21G x 1"	Posterior Capsule Polisher	10
M8604A	23G x 7/8"	Posterior Capsule Polisher	10
M8604B	21G x 7/8"	Posterior Capsule Polisher	10
M8605	27G x 7/8"	Silicone Tip Polisher/Manipulator	10
M8605A	23G x 7/8"	Silicone Tip Polisher/Manipulator	10
M8608	27G x 7/8"	Flexible Tipped Polisher	10
M8620	23G x 1"	Olive Tipped Cannula	10
M8620C	23G x 1"	Olive Tipped Cannula (Curved)	10
M8621	23G x 1"	Olive Tipped Cannula with Micro-etched Olive	10
M8622 NEW	25G x 7/8"	Olive Tipped Cannula (Angled) with Micro-etched Olive	10
M8645	23G x 7/8"	Posterior Scraper	10
M9272	23G x 1/2"	Lachrymal Cannula (Short Straight)	10
M9274	21G x 1 1/2"	Lachrymal Cannula (Curved)	10
M9275	26G x 3/4"	Lachrymal Cannula (Straight)	10
M9276	26G x 1 1/4"	Lachrymal Cannula (Curved)	10
M9277	25G x 1 1/8"	Lachrymal Cannula (Curved)	10
M9278	25G x 7/8"	Lachrymal Cannula (Straight)	10
M9279	19G/23G	Lachrymal Cannula / Dilator	10
M9279C	19G/23G	Lachrymal Cannula (Curved)	10
M9280	19G x 1"	Washout Cannula (Bishop Harmon)	10
M9281	16G x 1"	Washout Cannula	10

COLOR	GAUGE
Yellow 	30G
Grey 	27G
Brown 	26G
Orange 	25G
Purple 	24G
Blue 	23G
Black 	22G
Green 	21G
Yellow 	20G
Cream 	19G
Pink 	18G
White 	16G

All hub colours comply with ISO6009

KIT PACKERS

All Sterimedix products are available to kit packers, packed in single blisters, non-sterile. Drapes are provided in bulk, non-sterile.

In addition to Sterimedix products we also supply:
 Instrument Wipes
 Eye Shields
 Drains/Wicks

For more information,
 Contact us:
contact@eye4vision.com



www.eye4vision.com